



C10MR USB&Wifi HD

Linear Ultrasound Probe

Technical Specifications

INTRODUCTION



Improve patient safety with
ultra-high definition ultrasound
facial aesthetics

Nail Ultrasound

Finger nails revealed by
high-frequency ultrasound(US) in
patients with psoriasis arthritis


Mapping & Injecting

Visualize facial arteries in danger zones
With high-definition ultrasound for safer
injections

Evaluating Fillers

Confirm and characterize existing facial
fillers using ultra-high-frequency
ultrasound

MORE SPECS

Model	C10MR
Picture	
Probe Type	High Frequency Linear Probe
Software Systems	Android. IOS. Windows
Probe Button	<i>-Power on/off -Freeze/unfreeze</i>
*Connection	Wi-Fi USB
Wi-Fi Type	802.11n/20MHz/ 2.4G/5G

Frequency	18/24Mhz
Scan Depth	1~25mm
Length	L18*W5mm
Elements	192
DR	40-110Db
Dimensions	155x65x25mm
Weight	155g
Waterproof Grade	IPX8
Display Mode	B B+B B+M Color PW PDI
*Dual&Full screen Display	●

Software Control		
B GN(Gain)		•
Color GN(Gain)		•
D(Depth)		•
ENH(Enhancement)		•
DR(Dynamic Range)		•
H(harmonic)		•
Focus		•
8 TGC		•
Mid-Line		•
An note		•
Page Left/Right Up/Down		•
*Body Mark		•
Dicom		•
Image storage		JPEG, MP4, Dicom
Clip Loop		Automatic / manual 100-200-500-1000 frame
Frame rate		18~24 F/S
Presets	Thyroid, SmallParts, Pediatrics, Vascular, Carotid, Breast, MSK, VesselFlow, Nerve	
Measures	B: Length, Distance, Area/Circumference, Trace, Depth, Angle, B+M: HeartRate, Time, Distance, LVID B+PW: Velocity, Heart Rate(2), S/D, Depth, Time.	

Measure Magnify	●
*Built in cooling coating	●
*Working While charging	●
Battery	Built-in High capacity Lition batteries 2200 mAh,
Power Consumption	3 W (Live) / 1 W (freeze)
Continuous hours of scanning	4h 72h Ultra-long standby
Recharge time	3~4h
*Charging	Type-C USB Wireless charging
Language	Chinese, Ukrainian, German, English, Russian, France, Italian, Spanish, Portuguese (Brazil)

(“●” means applicable, “○” means inapplicable)



Search and download

“MY USG”

from App Store or Google Play
on your smartphones/tablets/laptops.
Contact us to get windows software.





Beijing Konted Medical Technology Co.,Ltd

Add: 1F, Building 3, No.27 Yongwang Road, Daxing Biological Pharmaceutical Industry Base, Daxing District,

102629, Beijing,China

Web:www.konted.cn

Tel: 8610-60219113

Tax: 8610-60219213

Service:info@konted.cn